

Neuron Activation Pod

BY IVAR GESTRANIUS & KEVIN LAHTINEN FOR LOOOK INDUSTRIES



Neuron Activation Pod

BY IVAR GESTRANIUS & KEVIN LAHTINEN FOR LOOOK INDUSTRIES

Peduce Stress & Anxiety

The **Neuron Activation Pod** is equipped with groundbreaking Finnish health technology with a singular focus on enhancing human well-being. This pod is a giant leap towards the office of the future. This is not just a silent capsule or a traditional sleep pod.

The Neuron Activation Pod transmits very low frequency vibrations in the body. The technology promotes the body's natural healing abilities, especially through the quality and power of sleep. The technology mechanically guides the user into a meditation-like state. The vibration directly affects the autonomic nervous system, so the mind does not have to make an effort to relax. The mechanisms of the nervous system are activated, and stress-related symptoms, such as anxiety, sleep disorders and, according to customer experiences, also pain are relieved. The preventive effect takes care of you every day.

Different effects are achieved in different power ranges. Vibration helps an athlete or an individual experiencing work-related stress to recover faster, and activation programmes prepare the body for future performance by increasing alertness.

The relaxation treatment is suitable for all ages and is safe to use. For instance, pacemakers, pain stimulators, and artificial metal joints do not prevent the use of Neuron Activation Pod.








Neuron Activation Pod Specifications

BY IVAR GESTRANIUS & KEVIN LAHTINEN FOR LOOOK INDUSTRIES

Description

Fully upholstered. Black castors as standard. Delivered with integrated technology, integrated touch screen and headphones. Product delivered with 60" cable. Headrest and footrest in black leather as standard. Removable and changeable cushions and headrest.

Please keep in mind that this product is delivered fully assembled. Due to its size, it may not be suitable for locations with narrow door openings, staircases or small elevators. However, it can be disassembled into three smaller modules. Remember to check the dimensions before ordering this product.



Specifications

Neuron Activation Pod with Privacy Screen

• Model: LI-NA-602S

Dimensions: 83"W x 32"D x 58"H

· Weight: 428 lbs

Frame:

- · Compression moulded birch
- Plywood
- · High quality Polyurethane foam
- Steel
- · Fabric

Seat:

- Plywood
- · PET felt
- · Fabric

Details:

- Steel
 - Aluminum

Castors:

- · Polyamide
- · Polyurethane
- · Stainless steel

Fabric:

Available with multicolor upholstery and vinyl. Please check our fabric list or price list for options.

Electrical:

Available with US plug and play technology.









Neuron Activation Pod Programs

BY IVAR GESTRANIUS & KEVIN LAHTINEN FOR LOOOK INDUSTRIES

You can use the Neuron Activation Pod's integrated touchscreen to easily select ready 10-40 minute programs for relaxation, activation or recovery. Take a powernap in the middle of a busy day or choose a relaxation program for sleep problems or stress. You can use activation programs to boost your metabolism, increase your alertness, or warm up your muscles. Recovery programs are ideal after sports or for reducing swelling, for example.

All programs promote sleep quality and contribute to both physical and mental recovery. The effects are individual, and regular use yields the best results.





Recovery

Recover after a workout



Relaxation

Relax after a stressful meeting




Powernap

Catch up on some Zzzzs




Activation

Activate yourself when you're feeling low on energy



Program Settings

- · Select your preferred duration
- · Choose music preference and volume
- Enable alarm



Adjustable Vibrations

- Fine tuning vibrations in 4 separate region
- · Adjusts specifically to each individual



Neuron Activation Pod Technology

BY IVAR GESTRANIUS & KEVIN LAHTINEN FOR LOOOK INDUSTRIES















